Implantologie | Implantology

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Innovative, in cooperation with Dr. Regina Becker – Düsseldorf newly developed and clinically tested special instruments for the maintenance therapy of implants as well as for the treatment of peri-implant inflammations (Mucositis or Peri-Implantitis). The KKD[®] *mf* “RB-LINE” instrument line has been especially designed and optimised for these scopes of application. The special raw material properties and the design of the instruments guarantee a careful and effective working.

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Big surface mirror focused on the clinical examination of the gingiva, the mouth floor- or pharynx area, the lingual tooth surfaces and prosthetic restorations. Its optical properties enable a photo documentation by the mirror.

Due to its elongated design the mirror KKD[®] *mf* RB-LINE 1 guarantees during the treatment (even without assistance) a safe retraction of e.g. the lingual soft tissues and a good optical control of the instrument application. The shank is furnished with a M 2.5 mm screw thread.

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The special zirconium nitride coating of the tips enables a very rich-in-contrast evaluation of the penetration depth of the probe in the peri-implant sulcus. Because of the higher vulnerability of peri-implant tissues to inflammations compared to the tooth, it's recommended to use always a clean probe for sulcus depth measurement. Spreading of plaque from proximate periodontal pockets should be absolutely avoided. The millimetre gradations allow an exact measurement of the sulcus depth. The shank is furnished with a M 2.5 mm screw thread.

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Universal instrument with zirconium nitride coated working tips made of titanium. The working tips allow a safe biofilm removal during the maintenance therapy on implants. The slight form and special angulation allow an atraumatic access even for thin and close implants. The special coating avoids a fast dulling of the tips, damages of the implant surface as well as metal abrasion. This instrument allows a more precise and more effective working compared to carbon instruments.

The KKD® *mf* RB-LINE TiZi 3 has been developed on the base of the well-known curette type Gracey and can be used for different implant diameters in all 4 quadrants. The instrument is only sharpened on the inner surface and rounded outside, thus avoiding an unintended thinning of the peri-implant soft tissue collar.

KKD® *mf* RB-LINE 4 Zirkonnitridbeschichteter Medikamentenspatel | Zirconium nitride coated spatula

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Zirconium nitride coated special spatula with round- and longitudinal holes for improved application of antimicrobials in the peri-implant sulcus. The slight design of the spatula allows the use of the instrument for all implant diameters.

1 St./ pc. REF 14475

HK RB-LINE 5 Titan-Reinigungsbürsten | Titanium Cleaning Brushes

300 - 1.200 min⁻¹

NANO
1 St./ pc. REF 14490

Pocket
1 St./ pc. REF 14491

Titanium cleaning brushes for innovative and precise biofilm removal on rough titanium surfaces.

The HK RB-LINE 5 brushes has been developed to remove precisely biofilms during peri-implantitis treatments even on hard accessible areas. Used in a more intensive way it's even possible to smooth rough implant surfaces.

Contrary to rotating diamond instruments with this brush a regular surface treatment without weakening of the implant is achieved. Therefore also a displacement of titanium particles into peri-implant tissue is minimised.