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- Stabilization of the film and protection against deformation
- Exact and reproducible positioning of the film and/or memory foil
- Only little assistance of the patient is necessary
- Distortion free and fully detailed radiographs
- Easy handling
- Significantly improved radiation protection due to RWT® window x-ray®
- Maximum radiation protection with RWT® SAFETY bite blocks

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Using the Right Angle-/Paralleling Technique, due to the decreased vertical angle, the apex is no longer overlapped by the zygomatic bone.

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Due to the perpendicular position of the film/image plate to the longitudinal axis of the object the dimensional distortions are reduced to a minimum. Anatomic accuracy of tooth size, length of canals etc. is assured. Geometrical lack of definitions is minimized.

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Cone cutting, as a result of axial misalignment of the long cone, manual film/image plate holding or swinging out of the x-ray unit, is eliminated.

Vergleich der RWT® Haltersysteme – made in Germany Comparison of the RWT® Holder Systems – made in Germany

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

The KKD® RWT® COLOR instrument sets together with a long cone (rectangular- or round type) constitute the implementation of the Right Angle-/Parallelizing Technique in x-ray diagnosis. The special stepped design of the newly developed KKD® RWT® COLOR aiming rings enables an easy and safe aiming of the long cone. The rings are universally suitable for all current x-ray units (e.g. Dentsply Sirona, Kavo, Dürr, Satelec, Planmeca etc.). Thus allowing that the x-ray beam is directed perpendicular to the recording plane of the chosen medium which has been positioned parallel to the long axis of the tooth.

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Whereas with rectangular long cones the exposed x-ray field is already reduced to the chosen film size, with round cones the complete surface corresponding to the cross-section of the cone is exposed. Using the RWT® window x-ray® COLOR instrument set the radiation protection of round and rectangular long cones is the same. The lead-lined aiming rings allow a limitation of the exposed x-ray field to the chosen film size and – depending on the chosen window – reduce the radiation exposure by up to 70 %. The scattered radiation intensity is reduced – the radiograph contrast is increased.

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

The KKD® RWT® COLOR instrument sets together with a long cone (rectangular- or round type) constitute the implementation of the Right Angle-/Paralleling Technique in conventional & digital x-ray diagnosis. The rings are universally suitable for all current x-ray units (e.g. Dentsply Sirona, Kavo, Dürr, Satelec, Planmeca etc.).

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
Precise aiming of all current round- and rectangular cones due to the special stepped design.

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Bite blocks at their best: sophisticated and more comfortable to handle.

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
Easy and fast assembling due to colour coding.

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
High-quality indicator arms made of stainless medical steel.

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
The new design of holes enables an even easier attaching and separating of the bite blocks to/ from the indicator arm.

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Bite blocks anterior and posterior with enhanced clamping part: image plates and films are gently hold in position.

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Blocks are more rounded, hence less edges. On the back appropriate sizes are indicated.

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
With the RWT® COLOR SAFETY bite blocks, the last remaining safety risk in the x-ray diagnosis can be eliminated. The radiation exposure, which is limited to the chosen film size due to the use of a rectangular long cone or RWT window x-ray® COLOR, emerges barely reduced behind the film and this in spite of the integrated lead foil. The Roentgenographic Directive demands of the responsible doctor or dentist to reduce the radiation exposure for the patient to a minimum respecting the necessary aims. Furthermore it demands: [...] parts of the body, which are not concerned, [...] have to be protected as far as possible against radiation exposure.

RWT® COLOR SAFETY bite blocks offer a radiation barrier directly behind the used roentgenographic medium and absorb efficiently with the integrated lead foil the emerging x-rays: the remaining outlet is limited to the less than 1 % of the applied dose.

Dosismessung bei Heliodent DS - 60 KV/ Measured Dose Of The Heliodent Ds - 60 Kv

Schaltzeit exposure [s]	Aufnahmemedium radiographic medium	applizierte Dosis applied dose [µGy]	Restdosis hinter Bißplatte remaining dose behind bite block [µGy]	Durchlassstrahlung outlet radiation [%]	Schwächung der Exposition attenuation of the exposure [%]
0,64	Ultraspeed	1625	12,70	0,78	99,22
0,32	Ektaspeed	858	6,27	0,73	99,27
0,25	Insight	636	4,88	0,77	99,23
0,10	Speicherfolie	251	1,90	0,77	99,23

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
Base: test report of the TUV South Germany dated 8.8.2003
 - measured values of the Heliodent DS -
 Measured values are comparable with other used radiation sources

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The Kentzler-Kaschner Dental company has used the new construction of the bite blocks to increase patient comfort, since the innovative soft border avoids and eases pressure pain on the soft tissues.

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Complete assortment for anterior-
posterior- and bitewing radiographs with
multiple bite blocks and indicator arms.
Sufficient to prepare and assemble
4 complete holders.

4 Führungsstangen, 4 Visierringe,
18 Bißplatten
4 indicator arms, 4 aiming rings,
18 bite blocks

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Special assortment for the use of the Right
Angle-/Paralleling Technique for bitewing
radiographs. Complete kit with multiple bite
blocks, ring and arm.

1 Führungsstange, 1 Visierring, 4 Bißplatten
1 indicator arm, 1 aiming ring, 4 bite blocks
REF 05431S

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Special assortment for the use of the Right
Angle-/Paralleling Technique for endodontic
radiographs. Complete kit with multiple bite
blocks, ring and arm.

1 Führungsstange, 1 Visierring, 4 Bißplatten
1 indicator arm, 1 aiming ring, 4 bite blocks
REF 05430S

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Abmessungen: ca. 600 x 150 mm

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

Space saving holder base for storing mounted RWT® holder – also suitable for other brands.

- Wall-mounting, incl. mounting parts
- stable design
- made in Germany
- wipe disinfectable.

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We recommend to use the self-adhesive one-way blocks XIOS of the company Dentsply Sirona. Just mount the block – put the sensor with protection sleeve onto the sticky side of the block and start.

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We also manufacture individual holder systems for sensors for the industry. Contact us and benefit from our Know-How!