


# CONTOUR™ POWDER LIQUID GOLD

OBERFLÄCHENINDIKATOR  
SURFACE INDICATOR



TAUB  
PRODUCTS



Zur optischen Kenntlichmachung der Oberfläche von Zähnen auf dem Gipsmodell und/oder Keramikrestauration, um Farbe, Form und Lichtbrechung des vorhandenen Zahnes auf den Zahnersatz naturgetreu zu imitieren.

- Schnelles Sichtbarmachen von Oberflächenungenauigkeiten
- Trocknet schnell und ist leicht wieder abzdampfen
- Optimale Möglichkeit um natürliche und schöne Restaurationen herzustellen

For optical indication of the teeth surface on plaster models and/or ceramic restorations in order to create a true-to-nature reproduction of form, colour and refraction of the existing tooth

- to see imperfections and irregularities very quickly
- dries quickly and can be easily removed using water / steam
- for highly natural and more beautiful restorations after glazing

Contour™ Powder Liquid Gold, 12 ml  
Verdünner/ thinner, 12 ml

REF 08575  
REF 08574

KENTZLER-KASCHNER DENTAL GmbH  
Geschäftsbereich TOPDENT  
Mühlgraben 36  
D-73479 Ellwangen/Jagst  
Telefon: +49 (7961) 90 73-0  
Fax: +49 (7961) - 5 20 31  
info@kkd-topdent.de

[www.kkd-topdent.de](http://www.kkd-topdent.de)




M070926 / 09.12.2021

# BRIGHT SPOT II

FLÜSSIGES OKKLUSIONS-  
PRÜFMITTEL (INDIKATOR)



LIQUID OCCLUSION INDICATOR



Zur einfachen Kenntlichmachung von Störkontakten bei schwer zu markierenden Oberflächen wie hochglanzpolierten Kroneninnenseiten und Metallflächen.

- Schichtstärke von nur ca. 7 micron
- spiegelt okklusale Störkontakte präzise wider
- wird auf die Gegen- oder benachbarte Seite aufgebracht und markiert schon beim geringsten Kontakt

Sollte nur auf trockenen Oberflächen verwendet werden. Nur auf versiegelten Gipsoberflächen verwenden um ein Eindringen der Flüssigkeit zu vermeiden und um eine einfachere Entfernung zu gewährleisten

For easy marking of high spots on hard accessible surfaces as high-glossed polished inner surfaces of crowns or metal surfaces.

- thin film of around 7 microns
- reflects occlusal high spots very precisely
- applied to the opposing or adjacent surface it transfers a mark at the lightest touch or contact

It is best when applied to dry surfaces. Use it only on sealed stone surfaces in order to avoid material from penetrating the surface and for easier removal of the film.

Bright Spot rot/ red, 15 ml  
REF **08656**

Verdünner/ thinner, 30 ml  
REF **08658**



[www.kkd-topdent.de](http://www.kkd-topdent.de)