


CONTOUR™ POWDER LIQUID GOLD

OBERFLÄCHENINDIKATOR
SURFACE INDICATOR

NEU
NEW

TP TAUB
PRODUCTS



Zur optischen Kenntlichmachung der Oberfläche von Zähnen auf dem Gipsmodell und/oder Keramikrestauration, um Farbe, Form und Lichtbrechung des vorhandenen Zahnes auf den Zahnersatz naturgetreu zu imitieren.

- Schnelles Sichtbarmachen von Oberflächenungenauigkeiten
- Trocknet schnell. Die flüssige Goldfarbe ermöglicht einen sehr dünnen Farbauftrag und trocknet dadurch schneller als herkömmlicher Puder oder Paste
- Optimale Möglichkeit um natürliche und schöne Keramikrestorationen herzustellen
- Leicht entfernbar/wieder abzdampfen

For optical indication of the teeth surface on plaster models and/or ceramic restorations in order to create a true-to-nature reproduction of form, colour and refraction of the existing tooth.

- to see imperfections and irregularities very quickly
- dries quickly. The liquid gold colour allows a thin layer, therefore it dries faster than conventional powder and paste
- for highly natural and more beautiful restorations after glazing
- can be easily removed using water / steam

Contour™ Powder, 12 ml
Verdünner/ thinner, 12 ml

REF **08575**
REF **08574**

www.kkd-topdent.de



M070924 / 9.12.2020